

Thermo Scientific Dharmacon Reverse Transfection Format (RTF): a rapid method for RNAi-based high-throughput studies of biological pathways

RNA interference (RNAi) screens are easier and more efficient with Dharmacon® RTF small interfering RNA (siRNA) libraries, which combine optimized broad-spectrum transfection reagents and potent Thermo Scientific Dharmacon SMARTpool siRNAs into simple, assay-ready and automation-compatible formats.

The availability of siRNA collections targeting related genes offers the possibility of high-throughput studies; however, managing large-scale transfection of siRNA reagents is very time consuming and costly. To facilitate high-throughput analyses, a strategy was developed to streamline the transfection workflow, permitting rapid and economical screening. Dharmacon RTF siRNA libraries are SMARTpool® siRNA reagents predispensed in a multiwell cell culture plate that is ready for resuspension and immediate use. The Dharmacon RTF siRNA libraries contain preselected groups of rationally designed SMARTpool siRNA reagents targeting genes confirmed to be relevant to a particular pathway or to be phylogenetically related to the indicated gene family as defined by the Gene Ontology Consortium (www.geneontology.org).

The distinguishing feature of Dharmacon RTF is that the SMARTpool reagents are provided in a pre-aliquoted format, such that cells are plated simultaneously into wells containing the individual SMARTpool silencing reagents rehydrated with a lipid-media mixture. Here we illustrate this novel transfection strategy in the functional analysis of targets implicated in clathrin-mediated endocytosis (CME) demonstrating Dharmacon RTF siRNA Libraries as invaluable tools for robust and reliable high-throughput screens.

Delivery: conventional transfection versus reverse transfection formats

Conventional transfection, also known as forward transfection (FT), is the most common technique for delivering siRNA into cells for gene silencing. This standard approach involves preplating cells ~1 d before treatment. On day 2, target-specific siRNA is complexed with a transfection reagent for delivery into the cells. Silencing is then assessed 24 to 48 h later (or more) depending on the target and the detection method. For multiple targets or high-throughput strategies, conventional FT requires significant handling of individual samples, such that screens with hundreds or thousands of genes quickly become laborious and cost-prohibitive.

To overcome this challenge, a new reverse transfection format (RTF) was developed. RTF differs from FT by virtue of streamlined transfection preparation and setup (Fig. 1). The SMARTpool reagents are pre-aliquoted as a Dharmacon RTF library in 96-well plates, supplied in quadruplicate and ready for a one-time transfection of 50 μ M siRNA into the cells of interest. Each plate also contains three negative and three positive control siRNAs and is shipped complete with DharmaFECT® Transfection Reagent and DharmaFECT Cell Culture Reagent (DCCR). Following a very simple protocol, the SMARTpool reagents are rehydrated in each well with a mixture of DCCR and DharmaFECT. After a short incubation period (30–90 min), cells are added to each well and incubated under standard conditions. Silencing is then assessed using an assay appropriate for the pathway of interest.




Figure 1: Workflow comparison of FT and RTF, illustrating the labor- and reagent-saving features. (a) FT siRNA-mediated silencing of specific CME genes affects internalization of radiolabeled Tfn. (b) RTF recapitulates phenotypes observed by FT of siRNAs targeting a collection of genes involved in CME. Fluorescent images of FITC-transferrin internalization in HeLa cells treated with target-specific siRNAs.

Under the conditions just described, RTF achieves efficient and specific target knockdown with minimal effect on cell viability, and at lower overall concentrations of transfection reagent and siRNA relative to FT. By eliminating complex and time-consuming manual processing, Dharmacon RTF reduces handling to merely rehydrating the pre-aliquoted SMARTpool siRNAs and adding the cells. Thus, Dharmacon RTF siRNA libraries minimize the potential for exposure to ribonucleases and provide means for reliable and robust targeted gene knockdown.

Application and experimental format

To demonstrate the usefulness of the RTF approach, it was compared to FT in a functional analysis of genes involved in the CME pathway. CME is an important process in higher eukaryotes for the internalization of nutrients, macromolecules, viruses and plasma membrane proteins from the extracellular environment¹ and is of interest to those studying host-viral pathogen interactions.

In previous work, 13 genes (CHC, β_2 -adaplin, dynamin II, CALM, Eps15, Eps15R, epsin, EEA1, Rab5a, Rab5b, Rab5c, CLCa and CLCb) were targeted in HeLa cells by SMARTpool reagents using conventional FT. These targets were chosen because of their purported involvement in CME and the availability of antibodies for detection². The consequence of RNAi-mediated gene knockdown was assessed by Western blot analysis, and—where antibodies were lacking—by monitoring expression levels of yellow fluorescent protein–tagged constructs. The effect of targeted gene knockdown was also assessed quantitatively by monitoring the internalization of radiolabeled transferrin (¹²⁵I Tfn) (Fig. 1a). Tfn is a ligand that, upon binding to its receptor, is endocytosed constitutively through clathrin-coated pits. In all cases, the corresponding SMARTpool reagent and at least one of the individual siRNAs reduced protein levels by >90% (all SMARTpool reagents and at least three of four individual siRNAs reduced protein levels by 75%)². It was determined that the analysis of an expanded set of 44 genes involved in CME (Table 1) represented a suitable test case for the feasibility of RTF-based delivery and silencing. The corresponding Dharmacon RTF siRNA collection was rehydrated, and then HeLa cells were plated at a density of 15,000 cells/well. Uptake of fluorescein isothiocyanate–labeled Tfn (FITC-Tfn) was assessed in a short-time-course assay at 37 °C to avoid the contribution of recycled ligand. The resulting phenotypes, determined by the presence or absence of labeled vesicles (endosomes) and labeled plasma membrane, were scored by visualization of FITC-Tfn internalization (Fig. 1b). For example, siRNA-mediated silencing of Eps15 had no effect on the endocytosis of FITCTfn, whereas knockdown of CHC and β_2 -adaplin resulted in strong visual phenotypes (Fig. 1b). These results were consistent with those of FT, demonstrating the

reproducibility of Dharmacon RTF siRNA libraries in rapid screens of large gene families.

Table 1

<i>AMPH</i>	<i>ARRB2</i>	<i>CLTC</i>	<i>HIP1</i>	<i>RAB11B</i>	<i>SYNJ1</i>
<i>AP1M1</i>	<i>CAV1</i>	<i>DNM2</i>	<i>HIP1R</i>	<i>RAB4A</i>	<i>SYNJ2</i>
<i>AP1M2</i>	<i>CAV2</i>	<i>DAB2</i>	<i>ITSN1</i>	<i>RAB4B</i>	<i>SYT1</i>
<i>AP2A1</i>	<i>CBL</i>	<i>EEA1</i>	<i>NEDD4</i>	<i>RAB5A</i>	<i>SYT2</i>
<i>AP2B1</i>	<i>CBLB</i>	<i>EPN1</i>	<i>NEDD4L</i>	<i>RAB5B</i>	
<i>AP2M1</i>	<i>CBLC</i>	<i>EPS15</i>	<i>NSF</i>	<i>RAB5C</i>	
<i>ARF6</i>	<i>CLTA</i>	<i>EPS15L1</i>	<i>PICALM</i>	<i>H3GLB1</i>	
<i>ARRB1</i>	<i>CLTB</i>	<i>GRB2</i>	<i>RAB11A</i>	<i>SH3GLB2</i>	

Table 1: Collection of SMARTpool reagents targeting 44 genes implicated in CME.

Conclusions

The application of RNAi-mediated silencing in high-throughput functional analyses can be a laborious endeavor by virtue of the number of samples to be processed. Dharmacon RTF siRNA Libraries represent an important advance in critical screening tools that combine optimized broad-spectrum transfection reagents and potent SMARTpool siRNAs into simple, assay-ready and automation-compatible formats. Depending on the screening assay, a screen of hundreds of targets could be completed within a matter of days.

As described here, Dharmacon RTF application distinguishes itself by permitting efficient transfection and rapid screens. The functional analysis of genes involved in CME was consistent with that achieved by FT, thus RTF provides a cost-effective strategy for quick, reliable high-throughput screening and target validation methods.

Additional Information

www.thermoscientific.com/dharmacon

References

1. Le Roy, C. & Wrana, J. Clathrin- and non-clathrin-mediated endocytic regulation of cell signalling. *Nat. Rev. Mol. Cell Biol.* 6, 112–126 (2005).
2. Huang, F. et al. Analysis of clathrin-mediated endocytosis of EGF receptor by RNA interference. *J. Biol. Chem.* 279, 16657–16661 (2004).

Contact Information (next page)

Troubleshooting

For technical information or troubleshooting contact
Thermo Scientific Genomics Tech Support:

In North America (US, Canada, Central/South America)
Techservice.genomics@thermofisher.com
+1 (800) 235-9880

In Europe (EU, Middle East, Africa)
Techservice.emea.genomics@thermofisher.com
(+44) 1372 840410

In Other Countries
www.thermo.com/dharmacondistributors

Literature Code: 00189-08-C-01-U

© 2010 Thermo Fisher Scientific Inc. All rights reserved.
alamarBlue is a trademark of Accumed International,
Inc. All other trademarks are the property of Thermo
Fisher Scientific Inc. and its subsidiaries.